

Certificate

To confirm the ability for continuos operation

of the MR series Magnetic Stirrers with hotplate

The MR series Magnetic Stirrers with hotplate feature overtemperature safety circuits according to DIN EN 61010-1:2001 and DIN EN 61010-2-010:2014 and therefore is designed for continuous operation.

This statement is made under the precondition that all units are operated in accordance with the operation manual and in accordance with good practice standards for safety in laboratories, rules for accident preventions, and compliance with directions on hazardous materials.

Schwabach, January 2019

Stefan Peters
Leitung Forschung & Entwicklung

Marcell Sarré
Leitung Qualitätsmanagement / Qualitätssicherung